Modello: PA-1

**Termometro a distanza ad
Infrarossi**

Manuale d'Uso

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Produttore:

Shenzhen Haoyuan Electronic Co., Ltd

Floor 3, workshop, Zhiyang Science Park, No. 014,

Qianzhangqi Road, Zhanggetang, Guanlan Street, Longhua

New District, Shenzhen, China

Model: PA-1

**Infrared Thermometer
User Manual**

Warning: It is the responsibility of the user to ensure the electromagnetic environment of the equipment so that the instrument can work normally.

Warning: Do not use this equipment near strong radiation sources, otherwise it may interfere with the normal operation of this instrument.

Warning: In addition to the transducers, accessories and cables sold by the manufacturer of this instrument or system as spare parts for internal components, the use of unspecified accessories, transducers and cables may result in increased emissions or immunity of the equipment or system Degree reduction.

Warning: This instrument should not be used close to or stacked with other equipment. If it must be used close to or stacked, it should be observed to verify that it can operate normally in its used configuration.

1. Infrared Thermometer

The infrared thermometer is mainly composed of an infrared temperature sensor, a cover, a display, a power supply circuit and a measurement circuit.

Application

The body temperature of the subject is displayed by measuring the heat radiation from the forehead.

Basic Working Principle

The infrared thermometer measures the infrared radiation exchange between the forehead of the human body to be measured and the appropriate correction value, and outputs the actual measurement data of the temperature of the human forehead. Finally, the actual measurement data of the measured object is displayed by measuring the heat radiation of the human forehead.

Product Features

- 1) High-precision infrared sensors, stable and reliable performance.
- 2) Measured values automatically saved.
- 3) Large-size LCD screen, high brightness backlight, clear and soft display.
- 4) Automatic shutdown saves energy and energy.
- 5) Clean and hygienic: non-contact forehead temperature measurement, measuring distance 1cm-15cm, no contact with human skin to prevent cross infection.
- 6) Simple and practical: One-button measurement, one-button memory check, easy operation.
- 7) Multiple sets of memory: 32 sets of measurement data can be stored for easy analysis and comparison.

Technical Indicators

Body temperature testing range	32.0°C~42.4°C
Measuring Accuracy	Temperature mode: $\pm 0.3^{\circ}\text{C}$.
Measuring distance	1 cm~15 cm
Operating environment temperature	Temperature range:10~ 40°C, humidity $\leq 85\%$
Transportation and storage temperature	Temperature range: -20°C ~+55°C humidity $\leq 93\%$
Power	DC9V
Auto power off	90 seconds
Battery Alert	Battery Low Alert
Backlight	High Brightness Backlight
Size	L 13.7 x W 5.3 x H 18.5 cm
Weight	160g (with battery)

2. Usage and Operation

Advice before Use

Infrared thermometers have high requirements on the surrounding environment, and correct use is the key to accurate testing.

- 1) This product is an infrared forehead thermometer that professionally measures the temperature of the human body. According to the difference of human skin, the measured temperature will be different.

- 2) When the room temperature changes greatly, the infrared frontal thermometer should be placed in the room for 20 minutes before use.
- 3) When the person to be measured comes from a place with a large difference in temperature from the measurement environment, they should stay in the test environment for at least 5 minutes, and wait for the temperature to be consistent with the environment before measuring.
- 4) The environment around the person to be tested must be stable. Do not test in places with large air flow such as fans and air outlets of air conditioners.
- 5) Do not use the instrument outdoors or in a place with strong sunlight.
If stored or used outside the temperature and humidity specifications specified above, the system may not meet the stated performance specifications.
- 6) Do not hold the front end of the thermometer during measurement.
- 7) The measurement site must not be covered by hair. If there is sweat, please wipe it with a dry towel and measure it later. Otherwise it will affect the measurement result.
- 8) Cold measurement of the forehead of a fever patient, sweating, and other cooling measures will make the measurement result low, and measurement should be avoided in this case.
- 9) Comparison of different temperature measurement methods. The body temperature measured by different temperature measurement methods will be different. The specific temperature differences are shown in the table below for reference.

Measurement Area	Normal Temperature
Anal temperature	36.6°C-38°C
Oral temperature	35.5°C-37.5°C
Underarm temperature	34.7°C-37.3°C
Cochlea temperature	35.8°C-38°C
Forehead temperature	35.8°C-37.8°C

Quick Installation

1. Install the battery.
2. The first time you use it, you need to put in a new battery, and then you can take a temperature measurement.
3. The front of the thermometer is aimed at the forehead, and the distance is 1 cm -15 cm. At the same time, make sure that there is no hair, sweat, cosmetics or hat covering, then press the measurement switch to start the measurement.
4. When the ambient temperature changes affect the forehead temperature or forehead sweating, please measure after measuring the earlobe. Also make sure that there is no hair, sweat, makeup or hat covering.

Note: The measured temperature value will vary depending on the person's skin color, thickness, and body part, and when the ambient temperature changes greatly, which is normal. The reason is that the more exposed human body parts are, the more affected by the ambient temperature.

Battery Replacement

The instrument uses a DC9V battery. Generally, it can be used continuously for more than 10000 times.

When a battery symbol appears on the top right of the screen and flashes, it means that the battery is insufficient, and the battery needs to be replaced as soon as possible.

- (1) Open the battery cover and replace the battery. Pay attention to the correct position of the positive and negative electrodes. The product is designed with a battery anti-reverse device. Even if the battery is installed in the wrong polarity, it will not cause damage to the product, but it will not work.
- (2) Do not use rechargeable batteries. Use disposable batteries. Use alkaline batteries whenever it is possible.
- (3) When the instrument is not working, the static current is extremely small. But when it is not used for a long time, it is recommended to remove the battery to prevent battery leakage from damaging the product.

Precautions

1. When opening the battery cover to replace the battery, pay special attention to the polarity of the battery. Misplacement may cause damage to the product.
2. When not in use for a long time, please remove the battery to prevent the thermometer from being damaged due to battery leakage.
3. Do not use if the battery is leaking or moldy.
4. Do not place the battery near a fire source or throw it into a fire to prevent the battery from exploding.
5. Do not store the battery in a high temperature or high humidity environment.
6. To avoid a short circuit, do not put the battery in the same


pocket with metal objects such as coins or keys or other containers that may shorten the battery life.

7. Do not dispose of used batteries arbitrarily, they should be classified and disposed according to local regulations.
8. Check the battery for power before using the thermometer every time.

Operating Instructions

Description of the thermometer:

1. Infrared sensor
2. Measurement port
3. LCD liquid crystal display
4. "+" key
5. "-" key
6. "M" key
7. Backlight keys
8. Buzzer button
9. Measurement button
10. Battery cover
11. Handle

Operation Steps

- (1) The thermometer is in the temperature measurement mode. Press and hold the power button to measure continuously.

If you don't need to measure for a long time, please release it after pressing the power button.

- (2) There is a counter in the lower left corner of the screen.
The count value is cycled from 1 to 32.
The counter will increase by 1 each time a valid measurement is performed, and then return to 1 after 32. If the measurement result is "LO" or "HI", it is called invalid measurement, the measurement result is not entered into the memory, and the counter is not incremented.
- (3) View historical data. After each test, the instrument will automatically record test data, which can record up to 32 historical measurement data.
If you want to view these data, press "M" key, then press "+" or "-", press "-" to display the latest measured temperature, and press "+" to record the farthest time. Press "+" or "-" again to cycle.

Key Description

- (1) Single-click the "M" key for historical data query, short press "+" to query the next set of temperature readings, and short press "-" to query the previous set of temperature readings.
A total of 32 sets of historical data can be viewed.
- (2) If the measurement result is "LO" or "HI", it is called invalid measurement.
The measurement result does not enter the memory and the counter does not increase.
- (3) Short press the backlight key to turn the backlight on and off.
- (4) Short press the speaker button to turn the sound on and off.

3. Symbols

Symbols	Description	Symbols	Description
	BF type equipment		Avoid raining
	Refers to the attached file		Reference to instructions

4. Warnings and Instructions

Read this manual carefully before use.

The ambient temperature of this product is 10°C~ 40°C, the optimal temperature is 25°C.

Do not use this product in an environment above 40°C or below 10°C.

Do not place this product too close to a charged object to avoid electric shock.

Do not use this product in an environment with relative humidity greater than 85%.

Do not place this product too close to the electromagnetic range (such as radio, mobile phone, etc.).

Do not expose this product to the sun or near the stove. Do not touch water.

Do not drop this product, do not use it if it is damaged.

Hair on the forehead, sweat, hat or scarf will affect the accuracy of the measurement data.

Please make sure that the measurement distance is not more than 15cm.

When forehead sweating or other reasons cause the forehead temperature to not normally respond to body temperature, please measure from the earlobe.

When cleaning is required, wipe the surface of the product lightly with alcohol.

In the event of a problem with the product, please contact your dealer and do not attempt to repair it yourself.

Make sure to remove the hair and sweat from your forehead before measuring. The use of this product is not a substitute for a doctor's diagnosis.

The protective glass outside the LCD frame is very important, and it is also the fragile part of the meter. Please use it with care. Do not charge non-rechargeable batteries and do not throw them into fire.

Do not expose the product to the sun, and do not touch the water.

5. Common Faults

If you encounter the following problems during use, please refer to the following guides to help solve the problem. If the problem is still not resolved, please call our after-sales service:

(1) The screen displays HI information.

In the "body temperature" mode, when the measured temperature exceeds 43 degrees, the product will display HI, first check whether the measured part of the human body has been in contact with an external heat source. If it is not confirmed, you need to call our after-sales service.

(2) L0 information is displayed on the screen.

In "body temperature" mode, when the measured temperature is lower than 32 degrees, the instrument will display L0. Follow the prompts in the form below to find the cause. If it is confirmed that the reason is not in the form, you need to call our after-sales service. The main reasons why the screen displays "L0".

The reason of L0 is displayed	Suggestion
Hair or sweat	Ensure there are no obstructions or sweat on the forehead.
Cold air blows the forehead	Ensure stable air during measurement
Forehead just cold compress	Measure 10 minutes after cold compress
The measuring distance is too far	The best measurement distance is 1 cm -15cm, and do not use the forehead gun against the skin

6. Maintenance

During the use of the product, please follow the prompts when you find the following conditions.

External dirt: Wipe the dirt with a clean soft cloth, or wipe with a cotton swab moistened with medical alcohol. Wiping with medical alcohol can also have sterilization and disinfection. Take care not to add too much water or alcohol to prevent damage to the inside of the product.

Internal dirt: The lens of the internal probe is an important device. Do not touch or press with your fingers or other objects, otherwise it will affect the preparation of the measured value. When the lens surface is found to be dirty, wipe the lens surface with a cotton swab moistened with 95% or more absolute alcohol.

Tip: Do not use 75% disinfectant alcohol to wipe the lens (there will be traces of residual water). Do not use other chemical wipes to wipe the lens (which will cause damage to the lens).

Storage: Keep in a dry and dark place out of direct sunlight. It is recommended to practice several times to become familiar with the measurement method and try not to change the product factory settings. The measurement results can't replace physician diagnosis.

During use if any problem, please contact the seller or manufacturer

Please dispose of the waste and residues at the end of the service life of the product according to local laws and regulations.

If there is a product failure, please do not disassemble and repair it by yourself, please contact the after-sales service phone. After the product reaches the end of its life, please send it to a suitable place to protect the environment. (Note: The disposal of wastes and residues should comply with the corresponding national laws and regulations.

7. Warranty

Carefully read all instructions before using the device.

This product is warranted to be free from defects in materials, workmanship and functionality for a period of 2 years starting from the date of purchase in case of receipt, for a period of 1 year starting from the date of purchase in case of invoice.

This warranty shall not apply to any products that have been subjected to misuse, negligence, accidental or abnormal conditions of operation, and does not covers the natural capacity fadeout.

Please refer to qualified and authorized people for the repair of the device. A non-professional repair may cause malfunction and risks and void your warranty rights.

8. After-sale service

Please contact your Retailer for any technical problem.

9. Transportation and Storage

It should be stored in a clean, cool, dry, and ventilated place. It should not be stored in the open air. It should not be stored in the same place with odorous and corrosive items.

Storage conditions: $-10^{\circ}\text{C} \sim 55^{\circ}\text{C}$, 30% RH \sim 93% RH, normal pressure.

CORRECT DISPOSAL OF THE BATTERY

According to the European Directive 2006/66/CE

This marking on battery indicates that the batteries in this product should not be disposed of with other household waste at the end of their working life. If batteries are not properly disposed of, these substances can cause harm to human

health or the environment. To protect natural resources and to promote material reuse, please separate batteries from other types of waste and recycle them through your local, free battery return system. When needed, please ask your local free battery return system to remove the battery from the case.

CORRECT DISPOSAL OF THE PRODUCT

According to the Directive 2012/19/UE

This marking on the product, accessories or literature indicates that the product and its electronic accessories should not be disposed of with other household waste at the end of their working life. To prevent possible harm to the environment of human health from uncontrolled waste disposal, please separated these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources. Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can take these items for environmentally safe recycling.

Imported and distributed by:

ASIA KINGSTON EU IND Ltd

3 Bratya Miladinovi Str, 1000 Sofia, Bulgaria

Manufacturer

Shenzhen Haoyuan Electronic Co., Ltd

Floor 3, workshop, Zhiyang Science Park, No. 014,

Qianzhangqi Road, Zhangetang, Guanlan Street, Longhua

New District, Shenzhen, China

MADE IN CHINA